

**A SYSTEM FOR CELL-BASED SCREENING**

(Case No. 97,022)

5

**BACKGROUND OF THE INVENTION**

**Field of The Invention**

This invention is in the field of introducing fluorescent reagents into or applying fluorescent reagents to cells and monitoring the fluorescence in the cells.

**Description of The Prior Art**

15

Drug discovery is a long, multiple step process involving the identification of specific disease targets, development of an assay based on a specific target, validation of the assay, optimization and automation of the assay to produce a screen, high throughput screening of compound libraries using the assay, hit validation and hit compound optimization. The output of this process is a lead compound that goes into preclinical and eventually clinical trials. In this process, the screening phase is distinct from the assay development phases and testing the efficacy of the compounds in living biological systems.

20  
25  
30  
35  
40  
45  
50  
55  
60  
65  
70  
75  
80  
85  
90  
95  
100  
105  
110  
115  
120  
125  
130  
135  
140  
145  
150  
155  
160  
165  
170  
175  
180  
185  
190  
195  
200  
205  
210  
215  
220  
225  
230  
235  
240  
245  
250  
255  
260  
265  
270  
275  
280  
285  
290  
295  
300  
305  
310  
315  
320  
325  
330  
335  
340  
345  
350  
355  
360  
365  
370  
375  
380  
385  
390  
395  
400  
405  
410  
415  
420  
425  
430  
435  
440  
445  
450  
455  
460  
465  
470  
475  
480  
485  
490  
495  
500  
505  
510  
515  
520  
525  
530  
535  
540  
545  
550  
555  
560  
565  
570  
575  
580  
585  
590  
595  
600  
605  
610  
615  
620  
625  
630  
635  
640  
645  
650  
655  
660  
665  
670  
675  
680  
685  
690  
695  
700  
705  
710  
715  
720  
725  
730  
735  
740  
745  
750  
755  
760  
765  
770  
775  
780  
785  
790  
795  
800  
805  
810  
815  
820  
825  
830  
835  
840  
845  
850  
855  
860  
865  
870  
875  
880  
885  
890  
895  
900  
905  
910  
915  
920  
925  
930  
935  
940  
945  
950  
955  
960  
965  
970  
975  
980  
985  
990  
995  
1000  
1005  
1010  
1015  
1020  
1025  
1030  
1035  
1040  
1045  
1050  
1055  
1060  
1065  
1070  
1075  
1080  
1085  
1090  
1095  
1100  
1105  
1110  
1115  
1120  
1125  
1130  
1135  
1140  
1145  
1150  
1155  
1160  
1165  
1170  
1175  
1180  
1185  
1190  
1195  
1200  
1205  
1210  
1215  
1220  
1225  
1230  
1235  
1240  
1245  
1250  
1255  
1260  
1265  
1270  
1275  
1280  
1285  
1290  
1295  
1300  
1305  
1310  
1315  
1320  
1325  
1330  
1335  
1340  
1345  
1350  
1355  
1360  
1365  
1370  
1375  
1380  
1385  
1390  
1395  
1400  
1405  
1410  
1415  
1420  
1425  
1430  
1435  
1440  
1445  
1450  
1455  
1460  
1465  
1470  
1475  
1480  
1485  
1490  
1495  
1500  
1505  
1510  
1515  
1520  
1525  
1530  
1535  
1540  
1545  
1550  
1555  
1560  
1565  
1570  
1575  
1580  
1585  
1590  
1595  
1600  
1605  
1610  
1615  
1620  
1625  
1630  
1635  
1640  
1645  
1650  
1655  
1660  
1665  
1670  
1675  
1680  
1685  
1690  
1695  
1700  
1705  
1710  
1715  
1720  
1725  
1730  
1735  
1740  
1745  
1750  
1755  
1760  
1765  
1770  
1775  
1780  
1785  
1790  
1795  
1800  
1805  
1810  
1815  
1820  
1825  
1830  
1835  
1840  
1845  
1850  
1855  
1860  
1865  
1870  
1875  
1880  
1885  
1890  
1895  
1900  
1905  
1910  
1915  
1920  
1925  
1930  
1935  
1940  
1945  
1950  
1955  
1960  
1965  
1970  
1975  
1980  
1985  
1990  
1995  
2000  
2005  
2010  
2015  
2020  
2025  
2030  
2035  
2040  
2045  
2050  
2055  
2060  
2065  
2070  
2075  
2080  
2085  
2090  
2095  
2100  
2105  
2110  
2115  
2120  
2125  
2130  
2135  
2140  
2145  
2150  
2155  
2160  
2165  
2170  
2175  
2180  
2185  
2190  
2195  
2200  
2205  
2210  
2215  
2220  
2225  
2230  
2235  
2240  
2245  
2250  
2255  
2260  
2265  
2270  
2275  
2280  
2285  
2290  
2295  
2300  
2305  
2310  
2315  
2320  
2325  
2330  
2335  
2340  
2345  
2350  
2355  
2360  
2365  
2370  
2375  
2380  
2385  
2390  
2395  
2400  
2405  
2410  
2415  
2420  
2425  
2430  
2435  
2440  
2445  
2450  
2455  
2460  
2465  
2470  
2475  
2480  
2485  
2490  
2495  
2500  
2505  
2510  
2515  
2520  
2525  
2530  
2535  
2540  
2545  
2550  
2555  
2560  
2565  
2570  
2575  
2580  
2585  
2590  
2595  
2600  
2605  
2610  
2615  
2620  
2625  
2630  
2635  
2640  
2645  
2650  
2655  
2660  
2665  
2670  
2675  
2680  
2685  
2690  
2695  
2700  
2705  
2710  
2715  
2720  
2725  
2730  
2735  
2740  
2745  
2750  
2755  
2760  
2765  
2770  
2775  
2780  
2785  
2790  
2795  
2800  
2805  
2810  
2815  
2820  
2825  
2830  
2835  
2840  
2845  
2850  
2855  
2860  
2865  
2870  
2875  
2880  
2885  
2890  
2895  
2900  
2905  
2910  
2915  
2920  
2925  
2930  
2935  
2940  
2945  
2950  
2955  
2960  
2965  
2970  
2975  
2980  
2985  
2990  
2995  
3000  
3005  
3010  
3015  
3020  
3025  
3030  
3035  
3040  
3045  
3050  
3055  
3060  
3065  
3070  
3075  
3080  
3085  
3090  
3095  
3100  
3105  
3110  
3115  
3120  
3125  
3130  
3135  
3140  
3145  
3150  
3155  
3160  
3165  
3170  
3175  
3180  
3185  
3190  
3195  
3200  
3205  
3210  
3215  
3220  
3225  
3230  
3235  
3240  
3245  
3250  
3255  
3260  
3265  
3270  
3275  
3280  
3285  
3290  
3295  
3300  
3305  
3310  
3315  
3320  
3325  
3330  
3335  
3340  
3345  
3350  
3355  
3360  
3365  
3370  
3375  
3380  
3385  
3390  
3395  
3400  
3405  
3410  
3415  
3420  
3425  
3430  
3435  
3440  
3445  
3450  
3455  
3460  
3465  
3470  
3475  
3480  
3485  
3490  
3495  
3500  
3505  
3510  
3515  
3520  
3525  
3530  
3535  
3540  
3545  
3550  
3555  
3560  
3565  
3570  
3575  
3580  
3585  
3590  
3595  
3600  
3605  
3610  
3615  
3620  
3625  
3630  
3635  
3640  
3645  
3650  
3655  
3660  
3665  
3670  
3675  
3680  
3685  
3690  
3695  
3700  
3705  
3710  
3715  
3720  
3725  
3730  
3735  
3740  
3745  
3750  
3755  
3760  
3765  
3770  
3775  
3780  
3785  
3790  
3795  
3800  
3805  
3810  
3815  
3820  
3825  
3830  
3835  
3840  
3845  
3850  
3855  
3860  
3865  
3870  
3875  
3880  
3885  
3890  
3895  
3900  
3905  
3910  
3915  
3920  
3925  
3930  
3935  
3940  
3945  
3950  
3955  
3960  
3965  
3970  
3975  
3980  
3985  
3990  
3995  
4000  
4005  
4010  
4015  
4020  
4025  
4030  
4035  
4040  
4045  
4050  
4055  
4060  
4065  
4070  
4075  
4080  
4085  
4090  
4095  
4100  
4105  
4110  
4115  
4120  
4125  
4130  
4135  
4140  
4145  
4150  
4155  
4160  
4165  
4170  
4175  
4180  
4185  
4190  
4195  
4200  
4205  
4210  
4215  
4220  
4225  
4230  
4235  
4240  
4245  
4250  
4255  
4260  
4265  
4270  
4275  
4280  
4285  
4290  
4295  
4300  
4305  
4310  
4315  
4320  
4325  
4330  
4335  
4340  
4345  
4350  
4355  
4360  
4365  
4370  
4375  
4380  
4385  
4390  
4395  
4400  
4405  
4410  
4415  
4420  
4425  
4430  
4435  
4440  
4445  
4450  
4455  
4460  
4465  
4470  
4475  
4480  
4485  
4490  
4495  
4500  
4505  
4510  
4515  
4520  
4525  
4530  
4535  
4540  
4545  
4550  
4555  
4560  
4565  
4570  
4575  
4580  
4585  
4590  
4595  
4600  
4605  
4610  
4615  
4620  
4625  
4630  
4635  
4640  
4645  
4650  
4655  
4660  
4665  
4670  
4675  
4680  
4685  
4690  
4695  
4700  
4705  
4710  
4715  
4720  
4725  
4730  
4735  
4740  
4745  
4750  
4755  
4760  
4765  
4770  
4775  
4780  
4785  
4790  
4795  
4800  
4805  
4810  
4815  
4820  
4825  
4830  
4835  
4840  
4845  
4850  
4855  
4860  
4865  
4870  
4875  
4880  
4885  
4890  
4895  
4900  
4905  
4910  
4915  
4920  
4925  
4930  
4935  
4940  
4945  
4950  
4955  
4960  
4965  
4970  
4975  
4980  
4985  
4990  
4995  
5000  
5005  
5010  
5015  
5020  
5025  
5030  
5035  
5040  
5045  
5050  
5055  
5060  
5065  
5070  
5075  
5080  
5085  
5090  
5095  
5100  
5105  
5110  
5115  
5120  
5125  
5130  
5135  
5140  
5145  
5150  
5155  
5160  
5165  
5170  
5175  
5180  
5185  
5190  
5195  
5200  
5205  
5210  
5215  
5220  
5225  
5230  
5235  
5240  
5245  
5250  
5255  
5260  
5265  
5270  
5275  
5280  
5285  
5290  
5295  
5300  
5305  
5310  
5315  
5320  
5325  
5330  
5335  
5340  
5345  
5350  
5355  
5360  
5365  
5370  
5375  
5380  
5385  
5390  
5395  
5400  
5405  
5410  
5415  
5420  
5425  
5430  
5435  
5440  
5445  
5450  
5455  
5460  
5465  
5470  
5475  
5480  
5485  
5490  
5495  
5500  
5505  
5510  
5515  
5520  
5525  
5530  
5535  
5540  
5545  
5550  
5555  
5560  
5565  
5570  
5575  
5580  
5585  
5590  
5595  
5600  
5605  
5610  
5615  
5620  
5625  
5630  
5635  
5640  
5645  
5650  
5655  
5660  
5665  
5670  
5675  
5680  
5685  
5690  
5695  
5700  
5705  
5710  
5715  
5720  
5725  
5730  
5735  
5740  
5745  
5750  
5755  
5760  
5765  
5770  
5775  
5780  
5785  
5790  
5795  
5800  
5805  
5810  
5815  
5820  
5825  
5830  
5835  
5840  
5845  
5850  
5855  
5860  
5865  
5870  
5875  
5880  
5885  
5890  
5895  
5900  
5905  
5910  
5915  
5920  
5925  
5930  
5935  
5940  
5945  
5950  
5955  
5960  
5965  
5970  
5975  
5980  
5985  
5990  
5995  
6000  
6005  
6010  
6015  
6020  
6025  
6030  
6035  
6040  
6045  
6050  
6055  
6060  
6065  
6070  
6075  
6080  
6085  
6090  
6095  
6100  
6105  
6110  
6115  
6120  
6125  
6130  
6135  
6140  
6145  
6150  
6155  
6160  
6165  
6170  
6175  
6180  
6185  
6190  
6195  
6200  
6205  
6210  
6215  
6220  
6225  
6230  
6235  
6240  
6245  
6250  
6255  
6260  
6265  
6270  
6275  
6280  
6285  
6290  
6295  
6300  
6305  
6310  
6315  
6320  
6325  
6330  
6335  
6340  
6345  
6350  
6355  
6360  
6365  
6370  
6375  
6380  
6385  
6390  
6395  
6400  
6405  
6410  
6415  
6420  
6425  
6430  
6435  
6440  
6445  
6450  
6455  
6460  
6465  
6470  
6475  
6480  
6485  
6490  
6495  
6500  
6505  
6510  
6515  
6520  
6525  
6530  
6535  
6540  
6545  
6550  
6555  
6560  
6565  
6570  
6575  
6580  
6585  
6590  
6595  
6600  
6605  
6610  
6615  
6620  
6625  
6630  
6635  
6640  
6645  
6650  
6655  
6660  
6665  
6670  
6675  
6680  
6685  
6690  
6695  
6700  
6705  
6710  
6715  
6720  
6725  
6730  
6735  
6740  
6745  
6750  
6755  
6760  
6765  
6770  
6775  
6780  
6785  
6790  
6795  
6800  
6805  
6810  
6815  
6820  
6825  
6830  
6835  
6840  
6845  
6850  
6855  
6860  
6865  
6870  
6875  
6880  
6885  
6890  
6895  
6900  
6905  
6910  
6915  
6920  
6925  
6930  
6935  
6940  
6945  
6950  
6955  
6960  
6965  
6970  
6975  
6980  
6985  
6990  
6995  
7000  
7005  
7010  
7015  
7020  
7025  
7030  
7035  
7040  
7045  
7050  
7055  
7060  
7065  
7070  
7075  
7080  
7085  
7090  
7095  
7100  
7105  
7110  
7115  
7120  
7125  
7130  
7135  
7140  
7145  
7150  
7155  
7160  
7165  
7170  
7175  
7180  
7185  
7190  
7195  
7200  
7205  
7210  
7215  
7220  
7225  
7230  
7235  
7240  
7245  
7250  
7255  
7260  
7265  
7270  
7275  
7280  
7285  
7290  
7295  
7300  
7305  
7310  
7315  
7320  
7325  
7330  
7335  
7340  
7345  
7350  
7355  
7360  
7365  
7370  
7375  
7380  
7385  
7390  
7395  
7400  
7405  
7410  
7415  
7420  
7425  
7430  
7435  
7440  
7445  
7450  
7455  
7460  
7465  
7470  
7475  
7480  
7485  
7490  
7495  
7500  
7505  
7510  
7515  
7520  
7525  
7530  
7535  
7540  
7545  
7550  
7555  
7560  
7565  
7570  
7575  
7580  
7585  
7590  
7595  
7600  
7605  
7610  
7615  
7620  
7625  
7630  
7635  
7640  
7645  
7650  
7655  
7660  
7665  
7670  
7675  
7680  
7685  
7690  
7695  
7700  
7705  
7710  
7715  
7720  
7725  
7730  
7735  
7740  
7745  
7750  
7755  
7760  
7765  
7770  
7775  
7780  
7785  
7790  
7795  
7800  
7805  
7810  
7815  
7820  
7825  
7830  
7835  
7840  
7845  
7850  
7855  
7860  
7865  
7870  
7875  
7880  
7885  
7890  
7895  
7900  
7905  
7910  
7915  
7920  
7925  
7930  
7935  
7940  
7945  
7950  
7955  
7960  
7965  
7970  
7975  
7980  
7985  
7990  
7995  
8000  
8005  
8010  
8015  
8020  
8025  
8030  
8035  
8040  
8045  
8050  
8055  
8060  
8065  
8070  
8075  
8080  
8085  
8090  
8095  
8100  
8105  
8110  
8115  
8120  
8125  
8130  
8135  
8140  
8145  
8150  
8155  
8160  
8165  
8170  
8175  
8180  
8185  
8190  
8195  
8200  
8205  
8210  
8215  
8220  
8225  
8230  
8235  
8240  
8245  
8250  
8255  
8260  
8265  
8270  
8275  
8280  
8285  
8290  
8295  
8300  
8305  
8310  
8315  
8320  
8325  
8330  
8335  
8340  
8345  
8350  
8355  
8360  
8365  
8370  
8375  
8380  
8385  
8390  
8395  
8400  
8405  
8410  
8415  
8420  
8425  
8430  
8435  
8440  
8445  
8450  
8455  
8460  
8465  
8470  
8475  
8480  
8485  
8490  
8495  
8500  
8505  
8510  
8515  
8520  
8525  
8530  
8535  
8540  
8545  
8550  
8555  
8560  
8565  
8570  
8575  
8580  
8585  
8590  
8595  
8600  
8605  
8610  
8615  
8620  
8625  
8630  
8635  
8640  
8645  
8650  
8655  
8660  
8665  
8670  
8675  
8680  
8685  
8690  
8695  
8700  
8705  
8710  
8715  
8720  
8725  
8730  
8735  
8740  
8745  
8750  
8755  
8760  
8765  
8770  
8775  
8780  
8785  
8790  
8795  
8800  
8805  
8810  
8815  
8820  
8825  
8830  
8835  
8840  
8845  
8850  
8855  
8860  
8865  
8870  
8875  
8880  
8885  
8890  
8895  
8900  
8905  
8910  
8915  
8920  
8925  
8930  
8935  
8940  
8945  
8950  
8955  
8960  
8965  
8970  
8975  
8980  
8985  
8990  
8995  
9000  
9005  
9010  
9015  
9020  
9025  
9030  
9035  
9040  
9045  
9050  
9055  
9060  
9065  
9070  
9075  
9080  
9085  
9090  
9095  
9100  
9105  
9110  
9115  
9120  
9125  
9130  
9135  
9140  
9145  
9150  
9155  
9160  
9165  
9170  
9175  
9180  
9185  
9190  
9195  
9200  
9205  
9210  
9215  
9220  
9225  
9230  
9235  
9240  
9245

Science Applications International Corporation (SAIC) 130 Fifth Avenue, Seattle, WA. 98109 describes an imaging plate reader. This system uses a CCD camera to image the whole area of a 96 well plate. The image is analyzed to calculate the total fluorescence per well for homogeneous assays.

5 Molecular Devices, Inc. describes a system (FLIPR) which uses low angle laser scanning illumination and a mask to selectively excite fluorescence within approximately 200 microns of the bottoms of the wells in standard 96 well plates in order to reduce background when imaging cell monolayers. This system uses a CCD camera to image the whole area of the plate bottom. Although this system measures signals originating from a cell monolayer at the bottom of the 10 well, the signal measured is averaged over the area of the well and is therefore still considered a homogeneous measurement, since it is an average response of a population of cells. The image is analyzed to calculate the total fluorescence per well for cell-based homogeneous assays.

15 Proffitt et. al. Cytometry 24: 204-213 (1996) describes a semiautomated fluorescence digital imaging system for quantifying relative cell numbers *in situ* in a variety of tissue culture plate formats, especially 96-well microtiter plates. The system consists of an epifluorescence inverted microscope with a motorized stage, video camera, image intensifier, and a microcomputer with a PC-Vision digitizer. Turbo Pascal software controls the stage and scans the plate taking multiple images per well. The software calculates total fluorescence per well, provides for daily calibration, and configures easily for a variety of tissue culture plate formats. 20 Thresholding of digital images and reagents which only fluoresce when taken up by living cells are used to reduce background fluorescence without removing excess fluorescent reagent. For example, in this system cells are pretreated with fluorescein diacetate (FDA) and loaded into 96-well plates.

25 Fluorescence microscopy of cells and tissues is well known in the art. A variety of methods have been developed to image fluorescent cells in a microscope and extract information about the spatial distribution and temporal changes occurring in these cells. An article by Taylor, et al. in American Scientist 80 (1992), p. 322-335 describes many of these methods and their applications. These methods have been designed and optimized for the preparation of a few specimens for high spatial and temporal resolution imaging measurements of distribution, 30 amount and biochemical environment of the fluorescent reporter molecules in the cells.

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25

Treating cells with dyes and fluorescent reagents and imaging the cells is well known in the art. There is also a considerable body of prior art related to genetic engineering of cells to produce fluorescent proteins, such as modified green fluorescent protein (GFP) as a reporter molecule. The green fluorescent protein (GFP) of the jellyfish *Aequorea victoria* is a protein with an excitation maximum at 395 nm and an emission maximum at 510 nm and does not require an exogenous factor. Uses of GFP for the study of gene expression and protein localization are discussed in more detail in papers by Chalfie et al. in Science 263 (1994), p. 12501-12504. Some properties of wild-type GFP are disclosed by Morise et al. in Biochemistry 13 (1974), p. 2656-2662, and Ward et al. in Photochem. Photobiol. 31 (1980), p. 611-615. An article by Rizzuto et al. in Curr. Biology 5 (1995), p. 635-642 discusses the use of wild-type GFP as a tool for visualizing subcellular organelles in cells. A paper by Kaether and Gerdes in Febs Letters 369 (1995), p. 267-271, reports the visualization of protein transport along the secretory pathway using wild-type GFP. The expression of GFP in plant cells is discussed by Hu and Cheng in Febs Letters 369 (1995), p. 331-334, while GFP expression in Drosophila embryos is described by Davis et al. in Dev. Biology 170 (1995), p. 726-729. U. S. Patent 5,491,084 describes expressing GFP from *Aequorea victoria* in cells as a reporter molecule fused to another protein of interest. PCT/DK 96/00052 relates to methods of detecting biologically active substances affecting intracellular processes by utilizing a GFP construct having a protein kinase activation site. Numerous references are related to GFP proteins in biological systems. For example, PCT/US94/10165 describes a system for isolating cells of interest utilizing the expression of a GFP like protein. PCT/GB96/00481 describes the expression of GFP in plants. PCT/US95/01425 describes modified GFP protein expressed in transformed organisms to detect mutagenesis. U. S. Patents 5,401,629 and 5,436,128 describe assays and compositions for detecting and evaluating the intracellular transduction of an extracellular signal. Recombinant cells that express cell surface receptors and contain reporter gene constructs that include transcriptional regulatory elements that are responsive to the activity of cell surface receptors are used.

## BRIEF DESCRIPTION OF THE DRAWINGS

- 5              Figure 1 shows a diagram of the components of the cell-based scanning system.
- 10             Figure 2 shows a schematic of the microscope subassembly
- Figure 3 shows the camera subassembly
- Figure 4 illustrates cell scanning system process
- Figure 5 illustrates a user interface showing major functions to guide the user
- Figure 6 illustrates data presentation on screen
- Figure 7 flow chart of processing step for the cell-based scanning system
- Figure 8 A-J illustrates the strategy of the Nuclear Translocation Assay
- Figure 9 is example data from a known inhibitor of translocation
- Figure 10 is example data from a known stimulator of translocation

00000000000000000000000000000000

**SUMMARY OF THE INVENTION**

The invention relates to a computer controlled optical-mechanical system for rapidly determining the distribution, environment, or activity of fluorescently labeled reporter molecules in cells for the purpose of screening large numbers of compounds for those that specifically affect particular biological functions. The invention involves:

- providing cells containing fluorescent reporter molecules in an array of locations,
- treating the cells in the array of locations with one or more reagents,
- imaging numerous cells in each location with a fluorescence microscope,
- converting the optical information into digital data,
- utilizing the digital data to determine the distribution, environment or activity of the fluorescently labeled reporter molecules in the cells and the distribution of the cells, and
- interpreting that information in terms of a positive, negative or null effect of the compound being tested on the biological function

The array of locations may be a microtiter plate or a microchip which is a microplate having cells in an array of locations. The invention includes an apparatus and a computerized method for acquiring data such as a digital frame grabber, processing, displaying and storing the data.

**DETAIL DESCRIPTION OF THE INVENTION**

Screening large numbers of compounds for activity with respect to a particular biological function requires preparing arrays of cells for parallel handling of cells and reagents. Standard 25 96 well microtiter plates which are 86 mm by 129 mm, with 6mm diameter wells on a 9mm pitch, are used for compatibility with current automated loading and robotic handling systems. The microplate is typically 20 mm by 30 mm, with cell locations that are 100-200 microns in dimension on a pitch of about 500 microns. Methods for making microplate are described in U.S. Serial No. 60/018,696 filed on May 30, 1996, assigned to the same assignee. This application is 30 incorporated herein by reference in its entirety. Microplates may consist of coplanar layers of materials to which cells adhere patterned with materials to which cells will not adhere, or etched

3-dimensional surfaces of similarly patterned materials. For the purpose of the following discussion, the terms 'well' and 'microwell' refer to a location in an array of any construction to which cells adhere and within which the cells are imaged. Microplates also include fluid delivery channels in the spaces between the wells. The smaller format of a microplate increases 5 the overall efficiency of the system by minimizing the quantities of the reagents, storage and handling during preparation and the overall movement required for the scanning operation. In addition, the whole area of the microplate can be imaged more efficiently, allowing a second mode of operation for the microplate reader as described later in this document.

Those skilled in this art will recognize a wide variety of ways to make and deliver 10 fluorescent reporter molecules to cells. For example, fluorescently labeled biomolecules such as proteins, phospholipids and DNA hybridizing probes, as well as fluorescent reagents specifically synthesized with particular chemical properties of binding or association have been used as fluorescent reporter molecules. Fluorescently labeled antibodies are particularly useful reporter molecules due to their high degree of specificity for attaching to a single molecular target in a mixture of molecules as complex as a cell, tissue or extract of either. 15

Fluorescently labeled reporter molecules are useful for determining the location, amount and chemical environment of the reporter. For example, whether the reporter is in a lipophilic membrane environment or in a more aqueous environment can be determined. The pH environment of the reporter can be determined. It can be determined whether a reporter having a chelating group is bound to an ion, such as  $\text{Ca}^{++}$ , or not. 20

Those skilled in the art will recognize a wide variety of ways to measure fluorescence. For example, some fluorescent reporter molecules exhibit a change in excitation or emission spectra, some exhibit resonance energy transfer where one fluorescent reporter loses fluorescence, while a second gains in fluorescence, some exhibit a loss (quenching) or 25 appearance of fluorescence, while some report rotational movements.

As mentioned earlier, see Description of the Prior Art above, a cell can be genetically engineered to express reporter molecules such as GFP coupled to a protein of interest. Using the method and systems of the present invention permits the determination of the site and extent of activity of the protein in the cells.

30 Figure 1 is a schematic diagram of the system for measuring the distribution, environment, or activity of fluorescent reporter molecules in cells. An inverted fluorescent

microscope 1 is a Zeiss Axiovert inverted fluorescence microscope which uses standard objectives with magnification of 1-100x to the camera, and a white light source (e.g. 100W mercury-arc lamp or 75W xenon lamp) with power supply 2. There is an XY stage 3 to move the plate 4 in the XY direction over the microscope objective. A Z-axis focus drive 5 moves the objective in the Z direction for focusing. A joystick 6 provides for manual movement of the stage in the XYZ direction. A high resolution digital camera 7 acquires images from each well or location on the plate. 8 is a camera power supply. 9 is an automation controller and 10 is a central processing unit. The PC 11 provides a display 12 and has associated software. The printer 13 provides for printing of a hard copy record. 14 are the microscope oculars.

Figure 2 is a schematic of the microscope assembly 1 showing in more detail the XY stage 3, Z-axis focus drive 5, joystick 6, light source 2, automation controller 9, oculars 14. 15 and 16 are cables to the computer and microscope, respectively. In addition, Figure 2 shows a 96 well microtiter plate 17 which is moved on the XY stage 3 in the XY direction. Light from the light source 2 passes through the PC controlled shutter 18 to a motorized filter wheel 19 with excitation filters 20. The light passes into filter cube 25 which has a dichroic mirror 26 and an emission filter 27. Excitation light reflects off the dichroic mirror to the wells in the microtiter plate 17 and fluorescent light 28 passes through the dichroic mirror 26 and the emission filter 27 and to the digital camera 7.

Figure 3 shows a schematic drawing of the camera assembly. The digital camera 7, which contains an automatic shutter for exposure control, receives fluorescent light 28 from the microscope assembly. A digital cable 30 transports digital signals to the computer. 31 is the camera power supply.

Figure 4 illustrates an alternative embodiment of the invention in which cells are in microwells 40 on a microplate 41. Typically the microplate is 20 mm by 30 mm as compared to a standard 96 well microtiter plate which is 86 mm by 129 mm. The microplate chamber 42 serves as a microfluidic delivery system for the addition of compounds to cells. The microplate 41 in the microplate chamber 42 is placed in an XY microplate reader 43. Digital data is processed as described above. The small size of this microplate system increases throughput, minimizes reagent volume and provides for the ability to control the distribution and placement of cells for fast and precise cell-based analysis. This information can be displayed on a PC screen 11 and made part of a bioinformatics data base 44. This data base is an important part of

the present invention because it not only permits storage and retrieval of data obtained through the methods of this invention, but also permits acquisition and storage of external data relating to cells. Figure 5 is a PC display which illustrates the operation of the software.

In addition to the advantages cited earlier, the higher density array of cells on a microplate allows the microplate reader to image the whole microplate at a low resolution of a few microns per pixel for high throughput and image particular locations on the microplate at a higher resolution of less than 0.5 microns per pixel for higher precision. These two resolution modes improve the overall throughput of the system.

Referring to Figure 7, at the start of an automated scan, the operator enters information 100 that describes the sample, specifies the filter settings and fluorescent channels to match the biological labels being used and the information sought and then adjusts the camera settings to match the sample brightness. For flexibility to handle a range of samples, the software next allows selection of various parameter settings used to identify nuclei, cytoplasm, different fluorescent reagents, cell selection settings and number of cells to be analyzed per well. These parameters are stored in the system's database for easy retrieval for each automated run. The system's interactive cell identification mode simplifies the selection of morphological 15 parameter limits such as the range of size, shape, and intensity of cells to be analyzed. The user specifies which wells of the plate the system will scan and how many fields or how many cells to analyze in each well. Depending on the setup mode selected by the user at step 101, the system either automatically pre-focuses the region of the plate to be scanned using an autofocus 20 procedure to "find focus" of the plate 102 or the user interactively pre-focuses 103 the scanning region by selecting three "tag" points which defines the rectangular area to be scanned by the system. A least-squares fit "focal plane model" is then calculated by the system from these tag points to estimate the focus of each well during a an automated scan. The focus of each well is 25 estimated by interpolating from the focal plane model during a scan.

During an automated scan, the software dynamically displays the status of scan in progress such as the number of cells that have been analyzed, the current well that is being analyzed, and images of each independent wavelength as they are acquired, and the result of the screen for each well as it is acquired. The plate 4 is scanned in a serpentine style as the software 30 automatically moves the motorized microscope XY stage 3 from well to well and field to field within each well of a 96-well plate. Those skilled in the programming art will recognize how to

adapt software for scanning of other microplate formats such as 24, 48, and 384 well plates. The scan pattern of the entire plate as well as the scan pattern of fields within each well are programmed. The system adjusts sample focus with an autofocus procedure 104 through the Z axis focus drive 5, controls filter selection via a motorized filter wheel 19 and acquires and 5 analyzes images of up to four different colors ("channels" or "wavelengths").

The autofocus procedure is called at a user selected frequency, typically for the first field in each well and then once every 4 to 5 fields within each well. The autofocus procedure calculates the starting Z-axis point by interpolating from the pre-calculated plane focal model.

Starting a programmable distance above or below this set point, the procedure moves the 10 mechanical Z-axis through a number of different positions, acquires an image at each, and finds the maximum of a calculated focus score that estimates the contrast of each image. The Z position of the image with the maximum focus score determines the best focus for a particular field. Those skilled in the art will recognize this as a variant of automatic focusing algorithms as described in the prior art in Harms et al. in Cytometry 5 (1984), p. 236-243, Groen et al. in 15 Cytometry 6 (1985), p. 81-91, and Firestone et al. in Cytometry 12 (1991), p. 195-206.

For acquisition of images, the camera's 7 exposure time is separately adjusted for each dye to ensure a high-quality image from each channel. Software procedures can be called, at the user's option, to correct for registration shifts between wavelengths by accounting for linear (X and Y) shifts between wavelengths before making any further measurements. The electronic shutter 18 is controlled so that sample photo-bleaching is kept to a minimum. Background shading and uneven illumination can also be corrected by the software using algorithms known 20 in the prior art.

In one channel, images are acquired of a primary marker 105 (typically cell nuclei counterstained with DAPI or PI fluorescent dyes) which are segmented ("identified") using an 25 adaptive thresholding procedure. The adaptive thresholding procedure 106 is used to dynamically select the threshold of an image for separating cells from the background. The staining of cells with fluorescent dyes can vary to an unknown degree across cells in a microtiter plate sample as well as within images of field of cells within each well of a microtiter plate. This variation can occur due to sample preparation and/or as a result of the nature of cell biology. A global threshold is calculated for the complete image to separate the cells from background. The global 30 adaptive techniques used in the system are variants of those described in prior art in Kittler et al.

in Computer Vision, Graphics, and Image Processing 30 (1985), p. 125-147, Ridler et al. in IEEE Trans. Systems, Man, and Cybernetics (1978), p. 630-632. These global adaptive thresholding techniques account for field to field variation. Another adaptive thresholding algorithm in the system utilizes local region thresholding in contrast to global image thresholding. Image analysis of local regions in the image leads to better overall segmentation since staining of cell nuclei (as well as other labeled components) can vary across an image. The global/local processing model increases throughout and improves cell segmentation. Using this global/local procedure, a reduced resolution image (reduced in size by a factor of 2 to 4) is first globally segmented (using adaptive thresholding) to find regions of interest in the image. These regions then serve as guides to more fully analyze the same regions at full resolution. A more localized threshold is then calculated (using again adaptive thresholding) for each region of interest.

The output of the segmentation procedure is a binary mask wherein the objects are white and the background is black. This binary image, also called a mask in the prior art, is used to determine if the field contains objects 107. The mask is labeled with a blob labeling algorithm whereby each object (or blob) has a unique number assigned to it. Morphological features, such as area and shape, of the blobs are used to differentiate blobs likely to be cells from those that are considered artifacts. The user pre-sets the morphological selection criteria by either typing in known cells morphological features or by using the interactive training utility. If objects of interest are found in the field, images are acquired for all other active channels 108, otherwise the stage is advanced to the next field 109 in the current well. Each object of interest is located in the image for further analysis 110. The software determines if the object meets the criteria for a valid cell nucleus 111 by measuring its morphological features (size and shape). For each valid cell, the XYZ stage location is recorded, a small image of the cell is stored, and features are measured 112. The cell scanning system can perform multiple tests on cellular samples by applying a number of analysis methods simultaneously to measure features at multiple wavelengths including:

1. the total fluorescent intensity within the cell nucleus for colors 1-4  
2. the area of the cell nucleus for color 1 (the primary marker)  
3. the shape of the cell nucleus for color 1 is described by three shape features:  
5           a) perimeter squared area  
              b) box area ratio  
              c) height width ratio  
4. the average fluorescent intensity within the cell nucleus for colors 1-4 (i.e. #1 divided  
by #2)  
5. the total fluorescent intensity of a ring outside the nucleus (see Figure 8) that  
represents fluorescence of the cell's cytoplasm (cytoplasmic mask) for colors 2-4  
10  
6. the area of the cytoplasmic mask  
7. the average fluorescent intensity of the cytoplasmic mask for colors 2-4 (i.e. #5  
divided by #6)  
8. the ratio of the average fluorescent intensity of the cytoplasmic mask to average  
fluorescent intensity within the cell nucleus for colors 2-4 (i.e. #7 divided by #4)  
9. the difference of the average fluorescent intensity of the cytoplasmic mask and the  
average fluorescent intensity within the cell nucleus for colors 2-4 (i.e. #7 minus #4)  
10. the number of fluorescent domains (also call spots, dots, or grains) within the cell  
nucleus for colors 2-4  
20

Features 1 through 4 are commonly used in a variety of image analysis applications and  
are well known in prior art. Features 5-9 have been developed specifically to provide  
measurements of a cell's fluorescent molecules within the local cytoplasmic region of the cell  
and the translocation (i.e. movement) of fluorescent molecules from the cytoplasm to the  
25 nucleus. These screen specific features are used for analyzing cells in microplates for the  
inhibition of nuclear translocation. Inhibition of nuclear translocation of transcription factors  
provides a novel approach to screening intact cells. An automated screen of an inhibitor of NF-  
 $\kappa$ B translocation has been successfully performed. A specific algorithm measures the amount of  
NF- $\kappa$ B probe in the nuclear region (feature 4) versus the local cytoplasmic region (feature 7) of  
30 each cell. Quantification of the difference between these two sub-cellular compartments  
provides a measure of cytoplasm-nuclear translocation (feature 9).

Feature 10 is used for counting of DNA or RNA probes within the nuclear region in colors 2-4. For example, DNA probes are commercially available for identifying the centromeres of specific chromosomes. Cells are three-dimensional in nature and when examined at a high magnification under a microscope one probe may be in-focus while another may be completely 5 out-of-focus. The cell screening system has a procedure for detecting three-dimensional probes in nuclei by acquiring images from multiple focal planes. The software moves the Z-axis motor drive 5 in small steps where the step distance is user selected to account for a wide range of different nuclear diameters. At each of the focal steps, an image is acquired. The maximum gray-level intensity from each pixel in each image is found and stored in a resulting maximum 10 projection image. The maximum projection image is then used to count the probes. The above algorithm work well in counting probes that are not stacked directly above or below another one. To account for probes stacked on top of each other in the Z-direction, users can select an option 15 to improve the counting of probes by analyzing probes in each of the focal planes acquired. In this mode, the scanning system performs the maximum plane projection algorithm as discussed above, detects probe regions of interest in this image, then further analyzes these regions in all the focal plane images.

After measuring cell features 112, the system checks if there are any unprocessed objects in the current field 113. If there are any unprocessed objects, it locates the next object 110 and checks if it meets the criteria for a valid cell nucleus 111, and measures its features. After it the 20 system has processed all the objects in the current field, it checks if it is done with the current plate 114. If the system is not done with the current plate, it check if it needs to find more cells in the current well 115. If it needs to find more cells in the current well it advances the XYZ stage to the next field within the current well 109 or it advances the stage to the next well 116 of the plate.

25 After a scan of a plate is complete, images and data can be reviewed with the system's image review, data review, and summary review facilities. All images, data, and settings from a scan are archived in the system's database for later review. Users can review the images alone of every cell analyzed by the system with an interactive image review procedure 117. The user can review data on a cell-by-cell basis using a combination of interactive graphs, a data spreadsheet 30 of features measured, and images of all the fluorescent channels of a cell of interest with the interactive cell-by-cell data review procedure 118. Graphical plotting capabilities are provided in

which data can be analyzed via interactive graphs such as histograms and scatter plots. Users can review summary data that are accumulated and summarized for all cells within each well of a plate with an interactive well-by-well data review procedure 119. Hard copies of graphs and images can be printed on a wide range of standard printers. All images and data are stored in a  
5 the system's database for archival and retrieval or for interface with a network laboratory management information system. Data can also be exported to other third-party statistical packages to tabulate results and generate other reports.

As a final phase of a complete scan, reports can be generated on one or more statistics of features measured. Users can generate a graphical report of data summarized on a well-by-well  
10 basis for the scanned region of the plate using an interactive report generation procedure 120. This report includes a summary of the statistics by well in tabular and graphical format and identification information on the sample. The report window allows the operator to enter comments about the scan for later retrieval. Multiple reports can be generated on many statistics and be printed with the touch of one button. Reports can be previewed for placement and data  
15 before being printed.

**EXAMPLE SCREEN 1 - Automated Screen for Compounds that Induce or Inhibit Nuclear Translocation of a DNA Transcription Factor**

Regulation of transcription of some genes involves activation of a transcription factor in the cytoplasm, resulting in that factor being transported into the nucleus where it can 5 initiate transcription of a particular gene or genes. This change in transcription factor distribution is the basis of a screen for the cell-based screening system to detect compounds which inhibit or induce transcription of a particular gene or group of genes. A general description of the screen is given followed by a specific example.

The distribution of the transcription factor is determined by labeling the nuclei with a 10 DNA specific fluorophore like Hoechst 33423 and the transcription factor with a specific fluorescent antibody. After autofocusing on the Hoechst labeled nuclei, an image of the nuclei is acquired in the cell-based screening system at 20x magnification and used to create a mask by one of several optional thresholding methods. The morphological descriptors of the regions defined by the mask are compared with the user defined parameters and valid nuclear masks are 15 identified and used with the following algorithm to extract transcription factor distributions. Each valid nuclear mask is eroded to define a slightly smaller nuclear region. The original nuclear mask is then dilated in two steps to define a ring shaped region around the nucleus, which represents a cytoplasmic region. The average antibody fluorescence in each of these two regions is determined, and the difference between these averages is defined as the NucCyt 20 Difference. Two examples of determining nuclear translocation are discussed below and illustrated in Figure 8A-J. Figure 8A illustrates an unstimulated cell with its nucleus 200 labeled with a blue fluorophore and a transcription factor in the cytoplasm 201 labeled with a green fluorophore. Figure 8B illustrates the nuclear mask 202 derived by the cell-based screening system. Figure 8C illustrates the cytoplasm 203 of the unstimulated cell imaged at a green 25 wavelength. Figure 8D illustrates the nuclear mask 202 is eroded (reduced) once to define a nuclear sampling region 204 with minimal cytoplasmic distribution. The nucleus boundary 202 is dilated (expanded) several times to form a ring that is 2-3 pixels wide that is used to define the cytoplasmic sampling region 205 for the same cell. Figure 8E further illustrates a side view which shows the nuclear sampling region 204 and the cytoplasmic sampling region 205. Using 30 these two sampling regions, data on nuclear translocation can be automatically analyzed by the

cell-based screening system on a cell by cell basis. Figure 8F-J illustrates the strategy for determining nuclear translocation in a stimulated cell. Figure 8F illustrates a stimulated cell with its nucleus 206 labeled with a blue fluorophore and a transcription factor in the cytoplasm 207 labeled with a green fluorophore. The nuclear mask 208 in Figure 8G is derived by the cell based screening system. Figure 8H illustrates the cytoplasm 209 of a stimulated cell imaged at a green wavelength. Figure 8I illustrates the nuclear sampling region 211 and cytoplasmic sampling region 212 of the stimulated cell. Figure 8J further illustrates a side view which shows the nuclear sampling region 211 and the cytoplasmic sampling region 212.

A specific application of this method has been used to validate this method as a screen.

10 A human chondrocyte cell line was plated in 96 well microtiter plates. Some rows of wells were titrated with IL-1 $\alpha$ , a known inducer of the nuclear transcription factor NF- $\kappa$ B. The cells were then fixed and stained by standard methods with a fluorescein labeled antibody to NF- $\kappa$ B, and Hoechst 33423. The cell-based screening system was used to acquire and analyze images from this plate and the NucCyt Difference was found to be strongly correlated with the amount of IL-1 $\alpha$  added to the wells as illustrated in Figure 9. In a second experiment, an antagonist to the receptor for IL-1, IL-1RA was titrated in the presence of IL-1 $\alpha$ , progressively inhibiting the translocation induced by IL-1 $\alpha$ . The NucCyt Difference was found to strongly correlate with this inhibition of translocation, as well as illustrated in Figure 10.

15 Additional experiments have shown that the NucCyt Difference gives consistent results over a wide range of cell densities and reagent concentrations, and can therefore be routinely used to screen compound libraries for specific nuclear translocation activity. Furthermore, the same method can be used with antibodies to other transcription factors, or GFP-transcription factor chimeras, in living and fixed cells, to screen for effects on the regulation of transcription of this and other genes.

20 Figure 6 is a representative display on a PC screen of data which was obtained in accordance with Example 1. Graph 1 300 plots the difference between the average antibody fluorescence in the nuclear sampling region and cytoplasmic sampling region, NucCyt Difference versus Well #. Graph 2 301 plots the average fluorescence of the antibody in the nuclear sampling region, NP1 average, versus the Well #. Graph 3 302 plots the average antibody fluorescence in the cytoplasmic sampling region, LIP1 average, versus Well #. The

software permits displaying data from each cell. For example, Figure 5 shows a screen display 406, the nuclear image 401, and the fluorescent antibody image 402 for cell #14.

NucCyt Difference referred to in graph 1 303 of Figure 6 is the difference between the average cytoplasmic probe (fluorescent reporter molecule) intensity and the average nuclear probe (fluorescent reporter molecule) intensity. The invention provides a computer means for converting the digital signal from the camera into this parameter and for plotting the parameter verses the well number.

NP1 average referred to in graph 2 304 of Figure 6 is the average of cytoplasmic probe (fluorescent reporter molecule) intensity within the nuclear sampling region. The invention provides a computer means for converting the digital signal from the camera into this parameter and for plotting the parameter verses the well number.

L1P1 average referred to in graph 3 305 of Figure 6 is the average probe (fluorescent reporter molecule) intensity within the cytoplasmic sampling region. The invention provides a computer means for converting the digital signal from the camera into this parameter and for plotting the parameter verses the well number.

DRAFT EDITION

**EXAMPLE SCREEN 2 - Automated Screen for Compounds that Induce or Inhibit Hypertrophy in Cardiac Myocytes**

Hypertrophy in cardiac myocytes has been associated with a cascade of alterations in gene expression and can be characterized in cell culture by an alteration in cell size, that is clearly visible in adherent cells growing on a coverslip. Preliminary experiments indicate that a screen can be implemented using the following strategy. Myocytes cultured in 96 well plates, can be treated with various compounds and then fixed and labeled with a fluorescent antibody to a cell surface marker and a DNA label like Hoechst. After focusing on the Hoechst labeled nuclei, two images are acquired, one of the Hoechst labeled nuclei and one of the fluorescent antibody. The nuclei are identified by thresholding to create a mask and then comparing the morphological descriptors of the mask with a set of user defined descriptor values. Local regions containing cells are defined around the nuclei. The limits of the cells in those regions are then defined by a local dynamic threshold operation on the same region in the fluorescent antibody image. A sequence of erosions and dilations is used to separate slightly touching cells and a second set of morphological descriptors is used to identify single cells. The area of the individual cells is tabulated in order to define the distribution of cell sizes for comparison with size data from normal and hypertrophic cells. In addition, a second fluorescent antibody to a particular cellular protein, such as one of the major muscle proteins actin or myosin can included.

Images of this antibody can be acquired and stored with the above images, for later review, to identify anomalies in the distribution of these proteins in hypertrophic cells, or algorithms can be developed to automatically analyze the distributions of the labeled proteins in these images.

### Additional Screens

Those skilled in the art will recognize a wide variety of distinct screens that can be developed. There is a large and growing list of known biochemical and molecular processes in cells that involve translocations or reorganizations of specific components within cells. The signaling pathway from the cell surface to target sites within the cell involves the translocation of plasma membrane-associated proteins to the cytoplasm. For example, it is known that one of the src family of protein tyrosine kinases, pp60c-src, translocates from the plasma membrane to the cytoplasm upon stimulation of fibroblasts with platelet-derived growth factor (PDGF). In contrast, some cytoplasmic components translocate from the cytoplasm to the plasma membrane upon stimulation of cells. For example, it is known that the GTP-binding proteins of the Rho family are maintained as cytoplasmic complexes with RhoGDI in resting cells, but are released and translocate to plasma membrane during cell activation. In addition, specific organelles, such as components of the cytoskeleton, nuclear envelope, chromatin, golgi apparatus, mitochondria, and endosomes are reorganized in response to specific stimuli. Finally, the targets for screening can themselves be converted into fluorescence-based reagents that report molecular changes including ligand-binding and post-translational modifications.

CONFIDENTIAL

References

- 1) C.S. Dulcey, J.H. Georger, V. Krauthamer, T.L. Fare, D.A. Stenger, and J.M. Calvert, "Deep UV Photochemistry of Chemisorbed Monolayers: Fabrication of Patterned Co- Planar Molecular Assemblies", *Science*, 1991, 252, 551.
- 5 2) J.M. Calvert, "Lithographically Patterned Self-Assembled Films", in **Thin Films, Vol. 20 - Organic Thin Films and Surfaces: Directions for the in the Nineties**, A. Ulman, Ed., Academic Press, Boston (1995), 109-141.
- 10 3) C.S. Dulcey, J.H. Georger, M.S. Chen, S. McElvany, C.E. O'Ferrall, V.I. Benezra, and J.M. Calvert, "Photochemistry and Patterning of Self-Assembled Monolayer Films Containing Aromatic Hydrocarbon Functional Groups", *Langmuir*, 1996, 12, 1638.
- 15 4) W.J. Dressick, C.S. Dulcey, M.S. Chen, and J.M. Calvert "Photochemical Studies of (Aminoethylaminomethyl) phenethyltrimethoxysilane (PEDA) Self-Assembled Monolayer Films", *Thin Solid Films*, 1996, in press.
- 5) D.A. Stenger, J.H. Georger, C.S. Dulcey, J.J. Hickman, A.S. Rudolph, T.B. Nielsen, S.M. McCort, and J.M. Calvert, "Coplanar Molecular Assemblies of Amino-and Perfluorinated Alkylsilanes: Characterization and Geometric Definition of Mammalian Cell Adhesion and Outgrowth", *Journal of the American Chemical Society*, 1992, 114, 8435.
- 20 6) B.J. Spargo, M.A. Testoff, T. Nielson, D.A. Stenger, J.J. Hickman, and A.S. Rudolph, "Spatially Controlled Adhesion, Spreading, and Differentiation of Endothelial Cells on Self-Assembled Molecular Monolayers, Proceedings of the National Academy of Sciences, 1994, 91, 11070.

- REF ID: A650
- 7) L.A. Chrisey, P.M. Roberts, V.I. Benezra, W.J. Dressick, C.S. Dulcey, and J.M. Calvert, "Selective Attachment of Synthetic DNA to Self-Assembled Monolayer Functionalized Surfaces", **Proceedings, Materials Research Society**, 1994, 330, 179.
- 5 8) L.A. Chrisey, C.E. O'Ferrall, B.J. Spargo, C.S. Dulcey, and J.M. Calvert "Fabrication of Patterned DNA Surfaces", **Nucleic Acids Research**, 1996, 24, 3040.
- 10 9) L.A. Chrisey, G.U. Lee, and C.E. O'Ferrall, "Covalent Attachment of Synthetic DNA to Self-Assembled Monolayer Films", **Nucleic Acids Research**, 1996, 24, 3031.
- 15 10) S.K. Bhatia, J.L. Teixeira, M. Anderson, L.C. Shriver-Lake, J.M. Calvert, J.H. Georger, J.J. Hickman, C.S. Dulcey, P.E. Schoen, and F.S. Ligler, "Fabrication of Surfaces Resistant to Protein Adsorption and Their Application to Two-Dimensional Protein Patterning", **Analytical Biochemistry**, 1993, 208, 197.
- 11) J.M. Calvert, G.S. Calabrese, J.F. Bohland, M.S. Chen, W.J. Dressick, C.S. Dulcey, J.H. Georger, Jr., J. Kosakowski, E.K. Pavelchek, K.W. Rhee, and L.M. Shirey, "Photoresist Channel-Constrained Electroless Metal Deposition on Ligating Self-Assembled Films", **Journal of Vacuum Science and Technology**, 1994, B12, 3884.
- 20 12) R. Kapur, B.J. Spargo, M.S. Chen, J.M. Calvert, and A.S. Rudolph, "Fabrication of 3-Dimensionally Textured Biomedical Polymers from Etched Silicon Substrates", **Journal of Biomedical Materials Research**, in press.

- 13) Harms, H. and Aus, H. M., "Comparison of Digital Focus Criteria for a TV Microscope System," **Cytometry** 5 (1984), p. 236-243
- 14) Groen, F. A., Young, I. T. , and Lighthart, G., "A Comparison of Different Focus Functions for Use in Autofocus Algorithms," **Cytometry** 6 (1985), p. 81-91
- 5 15) Firestone, L., Cook, K., Talsania, N., Preston, K., "Comparison of Autofocus Methods for Automated microscopy," **Cytometry** 12 (1991), p. 195-206
- 16) Kittler, J. and Illingworth, J., "Threshold Selection Based on a Simple Image Statistic," **Computer Vision, Graphics, and Image Processing** 30 (1985), p. 125-147,
- 10 17) Ridler, T. and Calvard, S. "Picture Thresholding Using an Iterative Selection Method," **IEEE Trans. Systems, Man, and Cybernetics** (1978), p. 630-632